



EasyTeleMed
HOME



THE APP FOR PATIENTS

EasyTeleMed Home is the **CE-certified mobile application of EasyTeleMed** (class II/A medical device according to EU MDR 2017/745) for **autonomous acquisition of vital parameters from the patient and automatic and immediate transmission to the certified EasyTeleMed EMR Plus** or to any other electronic medical record.

EasyTeleMed Home has a **user-friendly graphical interface**: through **audio-visual reminders**, it reminds the patient when to acquire measurements according to the personalized care plan and guides the patient in the correct execution of activities, providing


immediate feedback on the outcome. It is also possible to record measurements outside the care plan and answer questionnaires.

Acquisition of vital parameters is performed via **Bluetooth biomedical sensors** and is automatic, eliminating errors due to manual transcription (manual entry is however guaranteed). Transmission of collected data to the cloud is also automatic.


EasyTeleMed Home is available as a stand-alone version (mobile application) and as a kit version, consisting of a tablet with the application installed and a set of Bluetooth biomedical sensors.

KEY FEATURES

- ✓ Simple and intuitive graphical user interface to guide the patient in the autonomous execution of measurements
- ✓ Audio-visual reminders for planned activities
- ✓ Large number of Bluetooth biomedical sensors integrated
- ✓ Any Bluetooth sensor can be integrated
- ✓ Acquisition and storage of scheduled and extra-agenda measurements from the supplied medical sensors
- ✓ Manual input of measurements
- ✓ Integration with any electronic medical record
- ✓ Forward of received measurements to EasyTeleMed EMR Plus or electronic medical record via secure HTTPS connection
- ✓ Possibility to view the history of the acquired measurements
- ✓ Televisits
- ✓ Completely wireless solution, also usable on the move
- ✓ Support of evaluation questionnaires
- ✓ Push notifications for immediate updates from EasyTeleMed EMR Plus or any electronic medical record












HOW IT WORKS



MEASURABLE PARAMETERS

Possibility of integrating new vital parameters

-  HEART RATE
-  BLOOD PRESSURE
-  GLICEMIA
-  QUESTIONNAIRES
-  BODY WEIGHT
-  SPO2
-  BODY TEMPERATURE
-  AUSCULTATION
-  PHYSICAL ACTIVITY

AVAILABLE DEVICES

Possibility of integrating new Bluetooth sensors

BLOOD PRESSURE	A&D UA-651BLE, A&D UA-767PBT-CI, TAIDOC TD-3128
SPO2	GIMA OXY10, MIR SPIROTEL II, NONIN 3230, TAIDOC TD-8255
GLICEMIA	ACCU-CHEK GUIDE
BODY WEIGHT	A&D UC-351PBT-CI, A&D UC-352BLE
TEMPERATURE	TAIDOC TD-1241
AUSCULTATION	MINTTI SMARTH0-D2
PHYSICAL ACTIVITY	TICWATCH C2

